

# ninoSAFE Cyto

*No limit –  
Superior Low Noise*

ninoSAFE Cyto – A very sophisticated 3-filter cabinet.  
Equipped with circular, easy to change H14-prefilters with exceptionally low pressure drop. This means low service costs, very low noise level as well as very low energy consumption and accordingly very low CO<sub>2</sub> footprint. The circular design makes the H14-prefilters very compact which enable an excellent sitting position with good leg space. ninoSAFE Cyto is manufactured and tested according to EN12469 as well as DIN12980. ninoSAFE Cyto offers the same standard features and customization opportunities as ninoSAFE Class II range.



*Circular H14 prefilters*

- Product protection**
- Operator protection**
- Environmental protection**
- All air double HEPA-filtered**

Contaminated air   Filtered air   Ambient air

