

# ninoSAFE Class II Pro

## Superior Operator and Product Protection

Intuitive controller with large double touch displays equipped with MODBUS interface for reading and controlling. Optional remote support via ethernet allowing technical support on distance wherever in the world the cabinet is installed.

ninoSAFE Pro can be delivered with a long list of options and can be customized to fulfil any requirement.

The control panel has integrated control of optional electrically elevated support which means that everything is controlled from one spot.

Excellent energy saving dimmable LED light gives up to 2000 Lux evenly and shadow free distributed over the whole working area. The back wall is powder coated in white colour for best lumina-tion and to reduce illumination.

Latest fan and filter technology together with excellent engineering of the body gives superior noise levels of all sizes – from 46 dB(A).



**Product protection**  
**Operator protection**  
**Environmental protection**

Contaminated air   Filtered air   Ambient air

