

ninoSAFE XL

Superior Inner Dimensions

EN 12469 and optionally DIN 12980

ninoSAFE XL Safety Cabinet has extraordinary inner dimensions for use where a normal MSC not has sufficient inside space. Interior height is 1000 mm and interior depth is 900 mm. The complete front can be opened to easily fit in large devices, e.g. a pipetting robot. ninoSAFE XL is manufactured according to EN 12469 and can optionally be delivered as a three-filter cabinet with H14 prefilters according to DIN12980.

ninoSAFE XL has large internal dimensions and therefore of course also large external dimensions. However, due to the smart design it is possible to finally assemble the unit on site in case of narrow doorways.

ninoSAFE XL can also be delivered as ninoSAFE Class II B2 XL cabinet with 100% exhaust or as ninoSAFE Sterile XL for product protection.



Optional H14 pre filters according to DIN12980