

ninoSAFE Class III

Superior Operator Protection – Glove box

A proper MSC Class III cabinet designed and manufactured according to EN 12469. Stainless steel chamber with large glove holes for best ergonomics. HEPA H14-filters on both supply and exhaust air with double exhaust HEPA as an option. Can be equipped with gas tight valves for both supply and exhaust. Easy access to all technical parts from the front. Optional transfer hatch on one or two sides makes it possible to transfer samples or products in a secure way. Can also be linked to a MSC Class II cabinet.

